

Napsin A

IN VITRO DIAGNOSTIC DATASHEET

INTENDED USE : IN VITRO DIAGNOSTIC USE

This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy. it is Rabbit Monoclonal.

DESCRIPTION : Napsin A is a human aspartic proteinase found primarily in type II pneumocytes, alveolar macrophages and in the renal proximal tubules. Napsin A is more sensitive and specific in detecting lung adenocarcinomas than TTF-1. Differentiation of NSCLC into histologic types is important with the emergence of new targeted therapies. Literature suggest that a panel consisting of Napsin A, TTF-1, p63 and CK 5/6 is useful in distinguishing lung adenocarcinoma from squamous cell carcinoma.

CATALOG NO : PL2121

PL2121-R7 7 ML RTU 70 TEST

PL2121-R1 1 ML RTU 10 TEST

STAINING PATTERN : Cytoplasmic

PL2121-1 1 ML 1/100 1000 TEST

PL2121-0,1 0,1 ML 1/100 100 TEST

POSITIVE CONTROL : Lung Adenocarcinoma

VOLUME : 7 ml Ready to Use (7 ml of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing stabilizing protein and 0.015mol/L sodium azide.)

HOST : Rabbit

CLONE : ZR206

ANTIBODY CONCENTRATION : Not known

SPECIES REACTIVITY : Human. Others not tested.

EPITOPE : Not determined

MICROBIOLOGICAL STATE : This product is not sterile.

PRETREATMENT : Staining of formalin-fixed tissue sections requires treating the tissue sections in boiling 10mM edta buffer, pH 8.0, for 10-20 minutes followed by cooling at room temperature for 20 min.

PRIMARY ANTIBODY INCUBATION TIME : 60 minutes at Room Temperature

STAINING TIPS : If the staining is too light, use lower dilution or longer time. If the staining is too strong, check pretreatment, use higher dilution or shorter time.

STORAGE AND STABILITY : This product contains sodium azide and is stable for 24 months when stored at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent is not stored as recommended, performance must be validated by the user.

TROUBLESHOOTING : Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).