

OCT-2

IN VITRO DIAGNOSTIC DATASHEET

INTENDED USE: IN VITRO DIAGNOSTIC USE

This Rabbit Monoclonal product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.

DESCRIPTION: Oct-2 is a transcription factor of the POU homeo-domain family that binds to the lg gene octamer sites, regulating Bcell-specific genes. The oct-2 expression can be used as a Bcell lineage and differentiation marker. Germinal center Bcells, mantle B-cells, monocytoid B-cells, and plasma cells show high-level expression of Oct-2. Additionally, mantle cell lymphoma, follicular lymphoma, marginal zone lymphoma, plasmacytoma, Burkitt lymphoma, diffuse large cell lymphoma, diffuse large B-cell lymphoma, Hodgkin lymphoma display increased expression of Oct-2. Several studies of Oct-2 expression have shown a low-level expression in pre-B, T-cell, myelomonocytic, and epithelial cell lines. In contrast, all mature B-cell lines have high levels of expression.

 CATALOG NO:
 PL2022
 PL2022-R7
 7 ML RTU 70 TEST

 PL2022-R1
 1 ML RTU 10 TEST

 STAINING PATTERN:
 Nuclear
 PL2022-1
 1 ML 1/100 1000 TEST

 PL2022-0,1
 0,1 ML 1/100 100 TEST

POSITIVE CONTROL: Tonsil or lymph node

VOLUME: 7 ml Ready to Use (7 ml of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing

stabilizing protein and 0.015mol/L sodium azide.)

HOST: Rabbit

CLONE: PL-8

ANTIBODY CONCENTRATION: Not known

SPECIES REACTIVITY: Human. Others not tested.

MICROBIOLOGICAL STATE: This product is not sterile.

PRETREATMENT: Staining of formalin-fixed tissue sections requires treating the tissue sections in boiling 10mM citrate buffer, pH 6.0, for 10-20 minutes followed by cooling at room temperature for 20 min.

PRIMARY ANTIBODY INCUBATION TIME: 60 minutes at Room Temperature

STAINING TIPS: If the staining is too light, use lower dilution or longer time. If the staining is

too strong, check pretreatment, use higher dilution or shorter time.

STORAGE AND STABILITY: This product contains sodium azide and is stable for 24 months when stored

at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent

is not stored as recommended, performance must be validated by the user.

TROUBLESHOOTING: Please contact PatoLab Technical Support by e-mail (patolab@patolab.com.tr).

