

## HPV

### IN VITRO DIAGNOSTIC DATASHEET

#### INTENDED USE : IN VITRO DIAGNOSTIC USE

This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.

#### DESCRIPTION :

Papillomaviruses are a diverse group of DNA-based viruses. More than 100 different human papillomavirus (HPV) types have been characterized. Some HPV types cause benign skin warts, or papillomas, for which the virus family is named.

**CATALOG NO :** PL1826

PL1826-R7 7 ML RTU 70 TEST

PL1826-R1 1 ML RTU 10 TEST

**STAINING PATTERN :** Nuclear

PL1826-1 1 ML 1/200 2000 TEST

PL1826-0,1 0,1 ML 1/200 200 TEST

**POSITIVE CONTROL :** HPV Infected Tissue

**VOLUME :** 7 ml Ready to Use ( 7 ml of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing stabilizing protein and 0.015mol/L sodium azide. )

**HOST :** Mouse

**CLONE :** K1H8

**ANTIBODY CONCENTRATION :** 200ug / ml

**SPECIES REACTIVITY :** Human. Others not tested.

**EPITOPE :** Not Known

**MICROBIOLOGICAL STATE :** This product is not sterile.

**PRETREATMENT :** Staining of formalin-fixed tissue sections requires treating the tissue sections in boiling 10mM citrate buffer, pH 6.0, for 10-20 minutes followed by cooling at room temperature for 20 min.

**PRIMARY ANTIBODY INCUBATION TIME :** 30 minutes at Room Temperature

**STAINING TIPS :** If the staining is too light, use lower dilution or longer time. If the staining is too strong, check pretreatment, use higher dilution or shorter time.

**STORAGE AND STABILITY :** This product contains sodium azide and is stable for 24 months when stored at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent is not stored as recommended, performance must be validated by the user.

**TROUBLESHOOTING :** Please contact Patolab Technical Support by e-mail ( patolab@patolab.com.tr ).